



## Quadtools Order Form

Client Name: \_\_\_\_\_ Therapist: \_\_\_\_\_

Trial:  Application Approval # \_\_\_\_\_ Quote:  Purchase

Funder: Enable - MOH  or ACC  Accessable - MOH  or ACC  Other  \_\_\_\_\_

Delivery Address: \_\_\_\_\_

\_\_\_\_\_  
\_\_\_\_\_

**Cripper** Quantity \_\_\_\_\_

Mini 11"  Shorty 17"  Medium 23"  Original 29"  36" Long  48" Long

Lightweight  or Heavy Duty  Left  or Right

Note: No lightweight or heavy duty preference available for 36" & 48"

**Sip-N-Grip Reacher** Quantity \_\_\_\_\_

Shorty 17"  Original 29"

Lightweight  or Heavy Duty  Left  or Right

**Santoku Knife 6" Blade** Quantity \_\_\_\_\_

Lightweight  or Heavy Duty  Left  or Right

**Garden Shears** Quantity \_\_\_\_\_

Left  or Right

Forearm Length \_\_\_\_\_

Forearm Circumference \_\_\_\_\_

Wrist Thickness \_\_\_\_\_

Hand Thickness \_\_\_\_\_

**Note:** Please use the 'how to measure' instructions (on reverse side of this page) and the 'sizing chart' (available on our website) for exact measurements required. All Quadtools are made to measure so these measurements are critical to ensure a safe and appropriate fit.

## How To Measure Your Arm

The Crippler is custom-made for each person.

Measuring your arm is a critical part of the process to ensure a proper fit.



### Measuring Forearm Length

Measure the distance from the inside of your elbow to your wrist with the elbow and wrist bent.



Please measure to the nearest 1/8th inch.



### Measuring Forearm Circumference

Measure around the thickest part of your forearm near the elbow.

Please measure to the nearest 1/8th inch



### Measuring Wrist and Hand Thickness

Download and print 'The Crippler Sizing Chart' from our website [www.c1south.co.nz](http://www.c1south.co.nz) to measure your wrist and hand thickness.