



Treat-Lite Pressure Injury Prevention Cushion

A simple, cost-effective high risk air cushion that requires no pump or electricity and is machine washable. Ideal for any individual who needs to be seated for long durations of time and is at risk of pressure.


The Treat-Lite Pressure Care Cushion by Dan Medica South is a 3" high cushion with four layers of air-filled micro filaments. This new technology allows users the flexibility of using the cushion without the need for pumps or electricity.

Treat-Lite has a unique cross checked weave providing unsurpassable strength unlike any other mesh fibre on the market.

Features and Benefits:


- An air cushion without a pump. No adjustments or maintenance required
- High risk cushion that has been pressure tested to achieve between 55/60 mmHg, competing with some of the more expensive cushions on the market
- Lightweight with exceptional air-flow to prevent heat build-up
- Anti-bacterial and anti-fungal
- Machine washable in temperatures up to 40°C
- 2 year warranty

*16x18 can be rotated to fit 18x16 wheelchairs




Code	Size W x D	Weight Limit
TLite 1616	16" x 16"	125kg
TLite 1618	16" x 18"	125kg
TLite 1618*	18" x 16"	125kg
TLite 1818	18" x 18"	125kg
TLite 2020	20" x 20"	125kg

Without Treat-Lite®



With Treat-Lite®





Treat-Eezi Pressure Injury Prevention Overlays

Using the four layers of fine polyester fibres, Treat-Eezi overlays naturally and gently move with the user creating an almost zero friction area for a silent, comfortable sleep or seating surface along with the assurance that pressure injuries are given the best possible chance to breathe and heal. Proven to assist in healing Stage 2, 3 and 4 open wounds within four weeks!


Treat-Eezi Mattress Overlay

Fits single beds and hospital beds.
195cm x 90cm. 159kg weight limit.



Treat-Eezi Small Overlay

Targeted skin protection without the need for a full overlay.
100cm x 80cm. 159kg weight limit.



Treat-Eezi Seat Pad

40cm x 50cm (16" x 20").



The Treat-Eezi Seat Pad can be used over most wheelchair cushions or seating surfaces including recliner chairs, armchairs, car seats or office chairs for pressure relief. Minimises friction and shear forces, and distributes the weight bearing load by conforming to natural body contours for increased blood flow.

- Requires no electricity, pumps or maintenance
- Premium protection with exceptional pressure mapping results
- Improves blood flow and minimises shearing
- Reduces users temperature
- Lightweight, breathable and comfortable
- Can be bleached and washed in temperatures up to 80°C
- 2 year warranty

C1 SOUTH